



ZCCHC4 Antibody, Biotin conjugated

Product Code	CSB-PA887984LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5U6
Immunogen	Recombinant Human Zinc finger CCHC domain-containing protein 4 protein (25-229AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger CCHC domain-containing protein 4, ZCCHC4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZCCHC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.