





ZDHHC2 Antibody

Product Code	CSB-PA883437LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIJ5
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC2 protein (230-367AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC, IF; Recommended dilution: IHC:1:20-1:200, IF:1:50-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC2 (EC 2.3.1.225) (Reduced expression associated with metastasis protein) (Ream) (Reduced expression in cancer protein) (Rec) (Zinc finger DHHC domain-containing protein 2) (DHHC-2) (Zinc finger protein 372), ZDHHC2, REAM REC ZNF372
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ZDHHC2

Image



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA883437LA01HU at dilution of 1:100



Immunofluorescent analysis of HepG2 cells using CSB-PA883437LA01HU at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.