



# ZDHHC2 Antibody, FITC conjugated

<b>Product Code</b>	CSB-PA883437LC01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q9UIJ5
<b>Immunogen</b>	Recombinant Human Palmitoyltransferase ZDHHC2 protein (230-367AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Form</b>	Liquid
<b>Conjugate</b>	FITC
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Palmitoyltransferase ZDHHC2 (EC 2.3.1.225) (Reduced expression associated with metastasis protein) (Ream) (Reduced expression in cancer protein) (Rec) (Zinc finger DHHC domain-containing protein 2) (DHHC-2) (Zinc finger protein 372), ZDHHC2, REAM REC ZNF372
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	ZDHHC2
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.