



ZDHHC2 Antibody, HRP conjugated

Product Code	CSB-PA883437LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UIJ5
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC2 protein (230-367AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC2 (EC 2.3.1.225) (Reduced expression associated with metastasis protein) (Ream) (Reduced expression in cancer protein) (Rec) (Zinc finger DHHC domain-containing protein 2) (DHHC-2) (Zinc finger protein 372), ZDHHC2, REAM REC ZNF372
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ZDHHC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.