



ZDHHC3 Antibody, Biotin conjugated

Product Code	CSB-PA882141LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NYG2
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC3 protein (94-171AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC3 (EC 2.3.1.225) (Protein DHHC1) (Zinc finger DHHC domain-containing protein 3) (DHHC-3), ZDHHC3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	ZDHHC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.