



ZDHHC4 Antibody, HRP conjugated

Product Code	CSB-PA873612LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPG8
Immunogen	Recombinant Human Probable palmitoyltransferase ZDHHC4 protein (277-344AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable palmitoyltransferase ZDHHC4 (EC 2.3.1.225) (Zinc finger DHHC domain-containing protein 4) (DHHC-4) (Zinc finger protein 374), ZDHHC4, ZNF374
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZDHHC4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.