



ZDHHC6 Antibody, FITC conjugated

Product Code	CSB-PA862029LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H6R6
Immunogen	Recombinant Human Palmitoyltransferase ZDHHC6 protein (227-409AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Palmitoyltransferase ZDHHC6 (EC 2.3.1.225) (Transmembrane protein H4) (Zinc finger DHHC domain-containing protein 6) (DHHC-6) (Zinc finger protein 376), ZDHHC6, ZNF376
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZDHHC6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.