



ZDHHC6 Antibody, HRP conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA862029LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9H6R6 |
| Immunogen | Recombinant Human Palmitoyltransferase ZDHHC6 protein (227-409AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Palmitoyltransferase ZDHHC6 (EC 2.3.1.225) (Transmembrane protein H4) (Zinc finger DHHC domain-containing protein 6) (DHHC-6) (Zinc finger protein 376), ZDHHC6, ZNF376 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Cell Biology |
| Target Names | ZDHHC6 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |