



ZFAND2A Antibody, Biotin conjugated

Product Code	CSB-PA847647LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N6M9
Immunogen	Recombinant Human AN1-type zinc finger protein 2A protein (1-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	AN1-type zinc finger protein 2A, ZFAND2A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ZFAND2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.