




ZFP64 Antibody

Product Code	CSB-PA865094LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPA5
Immunogen	Recombinant Human Zinc finger protein 64 homolog, isoforms 1 and 2 protein (70-342AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 64 homolog, isoforms 1 and 2 (Zfp-64) (Zinc finger protein 338), ZFP64, ZNF338
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZFP64

Image



Western Blot

Positive WB detected in: Mouse liver tissue,
Mouse kidney tissue

All lanes: ZFP64 antibody at 4µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 75, 69, 73, 49, 47 kDa

Observed band size: 75 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.