



ZHX1-C8orf76 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA856963LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q96EF9 |
| Immunogen | Recombinant Human Zinc fingers and homeoboxes protein 1, isoform 2 protein (1-292AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc fingers and homeoboxes protein 1, isoform 2 (ZHX1-C8orf76 readthrough transcript protein), ZHX1-C8orf76, C8orf76 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Others |
| Target Names | ZHX1-C8orf76 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |