



ZHX1-C8orf76 Antibody, HRP conjugated

Product Code	CSB-PA856963LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EF9
Immunogen	Recombinant Human Zinc fingers and homeoboxes protein 1, isoform 2 protein (1-292AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc fingers and homeoboxes protein 1, isoform 2 (ZHX1-C8orf76 readthrough transcript protein), ZHX1-C8orf76, C8orf76
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZHX1-C8orf76
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.