



# ZKSCAN4 Antibody, Biotin conjugated

|                     |                                                                                                                                         |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA836168ED01HU                                                                                                                      |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                           |
| Uniprot No.         | Q969J2                                                                                                                                  |
| Immunogen           | Recombinant Human Zinc finger protein with KRAB and SCAN domains 4 protein (264-545AA)                                                  |
| Raised In           | Rabbit                                                                                                                                  |
| Species Reactivity  | Human                                                                                                                                   |
| Tested Applications | ELISA                                                                                                                                   |
| Form                | Liquid                                                                                                                                  |
| Conjugate           | Biotin                                                                                                                                  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, PH 7.4                                                        |
| Purification Method | >95%, Protein G purified                                                                                                                |
| Isotype             | IgG                                                                                                                                     |
| Clonality           | Polyclonal                                                                                                                              |
| Alias               | Zinc finger protein with KRAB and SCAN domains 4 (P373c6.1) (Zinc finger protein 307) (Zinc finger protein 427), ZKSCAN4, ZNF307 ZNF427 |
| Immunogen Species   | Homo sapiens (Human)                                                                                                                    |
| Research Area       | Epigenetics and Nuclear Signaling                                                                                                       |
| Target Names        | ZKSCAN4                                                                                                                                 |