



ZKSCAN7 Antibody, HRP conjugated

Product Code	CSB-PA026552LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0L1
Immunogen	Recombinant Human Zinc finger protein with KRAB and SCAN domains 7 protein (1-276AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein with KRAB and SCAN domains 7 (Zinc finger protein 167) (Zinc finger protein 448) (Zinc finger protein 64), ZKSCAN7, ZNF167 ZNF448 ZNF64
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZKSCAN7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.