



ZKSCAN8 Antibody

| | |
|----------------------------|---|
| Product Code | CSB-PA026568LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q15776 |
| Immunogen | Recombinant Human Zinc finger protein with KRAB and SCAN domains 8 protein (1-316AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc finger protein with KRAB and SCAN domains 8 (LD5-1) (Zinc finger protein 192), ZKSCAN8, ZNF192 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | ZKSCAN8 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |