



ZKSCAN8 Antibody, FITC conjugated

Product Code	CSB-PA026568LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15776
Immunogen	Recombinant Human Zinc finger protein with KRAB and SCAN domains 8 protein (1-316AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein with KRAB and SCAN domains 8 (LD5-1) (Zinc finger protein 192), ZKSCAN8, ZNF192
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZKSCAN8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.