



ZMYND19 Antibody, Biotin conjugated

Product Code	CSB-PA836202LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96E35
Immunogen	Recombinant Human Zinc finger MYND domain-containing protein 19 protein (1-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger MYND domain-containing protein 19 (Melanin-concentrating hormone receptor 1-interacting zinc finger protein) (MCH-R1-interacting zinc finger protein), ZMYND19, MIZIP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZMYND19