



ZMYND19 Antibody, Biotin conjugated

Product Code	CSB-PA836202LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96E35
Immunogen	Recombinant Human Zinc finger MYND domain-containing protein 19 protein (1-227AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger MYND domain-containing protein 19 (Melanin-concentrating hormone receptor 1-interacting zinc finger protein) (MCH-R1-interacting zinc finger protein), ZMYND19, MIZIP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZMYND19
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.