



ZNF19 Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA026567LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P17023 |
| Immunogen | Recombinant Human Zinc finger protein 19 protein (1-200AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc finger protein 19 (Zinc finger protein KOX12), ZNF19, KOX12 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | ZNF19 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |