




ZNF397 Antibody

Product Code	CSB-PA026720LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NF99
Immunogen	Recombinant Human Zinc finger protein 397 protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 397 (Zinc finger and SCAN domain-containing protein 15) (Zinc finger protein 47), ZNF397, ZNF47 ZSCAN15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZNF397

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA026720LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.