



ZNF610 Antibody, HRP conjugated

Product Code	CSB-PA026890LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N9Z0
Immunogen	Recombinant Human Zinc finger protein 610 protein (1-200AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 610, ZNF610
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF610
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.