




ZNF625 Antibody

Product Code	CSB-PA026903LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96I27
Immunogen	Recombinant Human Zinc finger protein 625 protein (1-306AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein 625, ZNF625
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF625

Image



Western blot

All lanes: ZNF625 antibody at 4µg/ml

Lane 1: Mouse brain tissue

Lane 2: Mouse lung tissue

Lane 3: Mouse ovarian tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 35, 43 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.