



ZNF738 Antibody, Biotin conjugated

Product Code	CSB-PA843334LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NE65
Immunogen	Recombinant Human Protein ZNF738 protein (1-137AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein ZNF738, ZNF738
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZNF738
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.