



ZPBP Antibody, Biotin conjugated

Product Code	CSB-PA866271LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BS86
Immunogen	Recombinant Human Zona pellucida-binding protein 1 protein (45-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zona pellucida-binding protein 1 (Inner acrosomal membrane IAM38) (Sp38), ZPBP, ZPBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZPBP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.