



ZBPB Antibody, FITC conjugated

Product Code	CSB-PA866271LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BS86
Immunogen	Recombinant Human Zona pellucida-binding protein 1 protein (45-351AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zona pellucida-binding protein 1 (Inner acrosomal membrane IAM38) (Sp38), ZBPB, ZBPB1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ZBPB