



ZSCAN16 Antibody, HRP conjugated

Product Code	CSB-PA862018LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4T2
Immunogen	Recombinant Human Zinc finger and SCAN domain-containing protein 16 protein (1-348AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and SCAN domain-containing protein 16 (Zinc finger protein 392) (Zinc finger protein 435), ZSCAN16, ZNF392 ZNF435
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZSCAN16