



ZSCAN16 Antibody, HRP conjugated

Product Code	CSB-PA862018LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H4T2
Immunogen	Recombinant Human Zinc finger and SCAN domain-containing protein 16 protein (1-348AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and SCAN domain-containing protein 16 (Zinc finger protein 392) (Zinc finger protein 435), ZSCAN16, ZNF392 ZNF435
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZSCAN16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.