



ZSCAN32 Antibody, Biotin conjugated

Product Code	CSB-PA026748LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NX65
Immunogen	Recombinant Human Zinc finger and SCAN domain-containing protein 32 protein (6-248AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and SCAN domain-containing protein 32 (Human cervical cancer suppressor gene 5 protein) (HCCS-5) (Zinc finger protein 434), ZSCAN32, ZNF434
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZSCAN32
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.