



ZSCAN9 Antibody, FITC conjugated

Product Code	CSB-PA026569LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15535
Immunogen	Recombinant Human Zinc finger and SCAN domain-containing protein 9 protein (1-394AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger and SCAN domain-containing protein 9 (Cell proliferation-inducing gene 12 protein) (PRD51) (Zinc finger protein 193), ZSCAN9, ZNF193
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZSCAN9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.