



ACOT11 Antibody

Product Code	CSB-PA855502LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WXI4
Immunogen	Recombinant Human Acyl-coenzyme A thioesterase 11 protein (35-254AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Acyl-coenzyme A thioesterase 11 (Acyl-CoA thioesterase 11) (EC 3.1.2.-) (Acyl-CoA thioester hydrolase 11) (Adipose-associated thioesterase) (Brown fat-inducible thioesterase) (BFIT), ACOT11, BFIT KIAA0707 THEA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	ACOT11
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.