



ACTR10 Antibody, HRP conjugated

Product Code	CSB-PA889156LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NZ32
Immunogen	Recombinant Human Actin-related protein 10 protein (165-417AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein 10 (Actin-related protein 11) (hARP11), ACTR10, ACTR11 ARP11
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ACTR10
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.