




ACTRT3 Antibody

Product Code	CSB-PA880135LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BYD9
Immunogen	Recombinant Human Actin-related protein T3 protein (1-372AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Actin-related protein T3 (ARP-T3) (Actin-related protein M1), ACTRT3, ARPM1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	ACTRT3

Image



Western blot

All lanes: ACTRT3 antibody at 2.5µg/ml

Lane 1: MCF-7 whole cell lysate

Lane 2: HCT116 whole cell lysate

Lane 3: K562 whole cell lysate

Lane 4: 293T whole cell lysate


Lane 5: HeLa whole cell lysate

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa



Immunohistochemistry of paraffin-embedded

human testis tissue using CSB-

PA880135LA01HU at dilution of 1:100

**Usage**

For Research Use Only. Not for use in diagnostic or therapeutic procedures.