




C11orf68 Antibody

Product Code	CSB-PA862008LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3H3
Immunogen	Recombinant Human UPF0696 protein C11orf68 protein (1-251AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UPF0696 protein C11orf68 (Basophilic leukemia-expressed protein Bles03) (Protein p5326), C11orf68, BLES03
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	C11orf68

Image



Western blot


All lanes: C11orf68 antibody at 3µg/ml + Mouse liver tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 28, 32 kDa

Observed band size: 28 kDa



Immunohistochemistry of paraffin-embedded human thyroid tissue using CSB-PA862008LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.