





C5orf22 Antibody

Product Code	CSB-PA675468LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q49AR2
Immunogen	Recombinant Human UPF0489 protein C5orf22 protein (101-420AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	UPF0489 protein C5orf22, C5orf22
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	C5orf22


Image



Western blot
All lanes: C5orf22 antibody at 8µg/ml + Mouse brain tissue
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 50, 21, 19 kDa
Observed band size: 50 kDa



Immunohistochemistry of paraffin-embedded human kidney tissue using CSB-PA675468LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA675468LA01HU at dilution of 1:100