



CCDC82 Antibody

Product Code	CSB-PA847633LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N4S0
Immunogen	Recombinant Human Coiled-coil domain-containing protein 82 protein (1-344AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Coiled-coil domain-containing protein 82, CCDC82
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CCDC82
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.