



CFHR2 Antibody, HRP conjugated

Product Code	CSB-PA13817B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36980
Immunogen	Recombinant Human Complement factor H-related protein 2 protein (20-269AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Complement factor H-related protein 2 (FHR-2) (DDESK59) (H factor-like 3) (H factor-like protein 2), CFHR2, CFHL2 FHR2 HFL3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CFHR2