




CHAD Antibody

Product Code	CSB-PA005310LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15335
Immunogen	Recombinant Human Chondroadherin protein (23-359AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chondroadherin (Cartilage leucine-rich protein), CHAD, SLRR4A
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	CHAD

Image



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA005310LA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.