



CXADR Antibody, FITC conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA006237LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P78310 |
| Immunogen | Recombinant Human Coxsackievirus and adenovirus receptor protein (259-365AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Coxsackievirus and adenovirus receptor (CAR) (hCAR) (CVB3-binding protein) (Coxsackievirus B-adenovirus receptor) (HCVADR), CXADR, CAR |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Microbiology |
| Target Names | CXADR |