





DDX31 Antibody

Product Code	CSB-PA872517LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H8H2
Immunogen	Recombinant Human Probable ATP-dependent RNA helicase DDX31 protein (560-851AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Probable ATP-dependent RNA helicase DDX31 (EC 3.6.4.13) (DEAD box protein 31) (Helicain), DDX31
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DDX31

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA872517LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human small intestine tissue using CSB-PA872517LA01HU at dilution of 1:100



Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.