



DPY30 Antibody, Biotin conjugated

Product Code	CSB-PA861193LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9C005
Immunogen	Recombinant Human Protein dpy-30 homolog protein (1-99AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein dpy-30 homolog (Dpy-30-like protein) (Dpy-30L), DPY30
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	DPY30
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.