



EDEM2 Antibody, Biotin conjugated

Product Code	CSB-PA861146LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV94
Immunogen	Recombinant Human ER degradation-enhancing alpha-mannosidase-like protein 2 protein (22-300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ER degradation-enhancing alpha-mannosidase-like protein 2, EDEM2, C20orf31 C20orf49
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	EDEM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.