



ELAVL1 Antibody, Biotin conjugated

Product Code	CSB-PA618993HD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15717
Immunogen	Recombinant Human ELAV-like protein 1 protein (2-326AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ELAV-like protein 1 (Hu-antigen R) (HuR), ELAVL1, HUR
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ELAVL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.