



EMILIN1 Antibody, Biotin conjugated

Product Code	CSB-PA896753LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6C2
Immunogen	Recombinant Human EMILIN-1 protein (1-347AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	EMILIN-1 (Elastin microfibril interface-located protein 1) (Elastin microfibril interfacier 1), EMILIN1, EMI
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	EMILIN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.