



ENY2 Antibody

Product Code	CSB-PA889067LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NPA8
Immunogen	Recombinant Human Transcription and mRNA export factor ENY2 protein (1-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcription and mRNA export factor ENY2 (Enhancer of yellow 2 transcription factor homolog), ENY2
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ENY2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.