



RTL8A Antibody

Product Code	CSB-PA871589LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BWD3
Immunogen	Recombinant Human Retrotransposon Gag-like protein 8A protein (1-113AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Retrotransposon Gag-like protein 8A (Mammalian retrotransposon derived protein 8A), RTL8A, CXX1B FAM127B MAR8A
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RTL8A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.