





FAM136A Antibody

Product Code	CSB-PA836190LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96C01
Immunogen	Recombinant Human Protein FAM136A protein (2-138AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:500-1:2000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM136A, FAM136A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM136A


Image



Western blot
All lanes: FAM136A antibody at 8µg/ml
Lane 1: Hela whole cell lysate
Lane 2: Caco-2 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 16 kDa
Observed band size: 16 kDa



Immunohistochemistry of paraffin-embedded human cervical cancer using CSB-PA836190LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA836190LA01HU at dilution of 1:100