



FAM71C Antibody, Biotin conjugated

Product Code	CSB-PA843337LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8NEG0
Immunogen	Recombinant Human Protein FAM71C protein (1-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein FAM71C, FAM71C
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM71C
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.