



CIAO2A Antibody, FITC conjugated

Product Code	CSB-PA872502LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H5X1
Immunogen	Recombinant Human Cytosolic iron-sulfur assembly component 2A protein (1-160AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytosolic iron-sulfur assembly component 2A, MIP18 family protein FAM96A, CIAO2A, CIA2A, FAM96A
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FAM96A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.