





FBXW10 Antibody

Product Code	CSB-PA726395LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5XX13
Immunogen	Recombinant Human F-box/WD repeat-containing protein 10 protein (701-1000AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	F-box/WD repeat-containing protein 10 (F-box and WD-40 domain-containing protein 10) (Ubiquitin ligase-specificity factor), FBXW10
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	FBXW10

Image



Immunohistochemistry of paraffin-embedded human colon cancer using CSB-PA726395LA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA726395LA01HU at dilution of 1:100

