



Fcgr3 Antibody, Biotin conjugated

Product Code	CSB-PA357412LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P08508
Immunogen	Recombinant Mouse Low affinity immunoglobulin gamma Fc region receptor III protein (31-215AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Low affinity immunoglobulin gamma Fc region receptor III (IgG Fc receptor III) (Fc-gamma RIII) (FcRIII) (CD antigen CD16), Fcgr3
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Fcgr3